


LB42708

Chemical Properties

CAS No. : 226929-39-1
 Formula: C₃₀H₂₇BrN₄O₂
 Molecular Weight: 555.47
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	LB42708 is an orally active farnesyltransferase (FTase) inhibitor (IC ₅₀ : 0.8/1.2/2.0 nM toward H/N/K-ras).
Targets(IC ₅₀)	Apoptosis,Transferase
In vitro	In mice treated with lipopolysaccharide (LPS), LB42708 (12.5 mg/kg, intraperitoneally) inhibited the production of nitric oxide (NO), prostaglandin E ₂ (PGE ₂), tumor necrosis factor-alpha (TNF-α), and interleukin-1 beta (IL-1β), and suppressed the progression of collagen-induced arthritis (CIA). Additionally, in tumors with either wild-type or mutated Ras, LB42708 (20 mg/kg/day, intraperitoneally) inhibited tumor growth and angiogenesis.
In vivo	LB42708 exerts an irreversible inhibitory effect on growth and induces apoptosis in rat intestinal epithelial cells altered by H-ras and K-ras. It inhibits VEGF-induced tumor angiogenesis in tumor-associated endothelial cells by blocking the Ras-dependent MAPK and PI3K/Akt signaling pathways. Additionally, in the murine macrophage lineage RAW264.7 cells, LB42708 significantly inhibits the process of lipopolysaccharide + IFN-γ induced intracellular farnesylation of the protein p21ras. In immune-activated osteoblasts and macrophages, LB42708 suppresses the production of nitric oxide synthase, cyclooxygenase-2, TNF-α, IL-1β, NO, and PGE(2). By inhibiting IKK activity, LB42708 also represses the activation of NF-κB and the activity of the iNOS promoter.
Kinase Assay	Biochemical Assessment of PDGFRα Kinase Activity: Chinese hamster ovary (CHO) cells are transiently transfected with mutated or wild type PDGFRα constructs and treated with various concentrations of Crenolanib. Experiments involving recombinant DNA are performed using biosafety level 2 conditions in accordance with guidelines. Protein lysates from cell lines are prepared and subjected to immunoprecipitation using anti-PDGFRα antibodies followed by sequential immunoblotting for PDGFRα. Densitometry is performed to quantify drug effect using Photoshop software, with the level of phosphor- PDGFRα normalized to total protein. Densitometry and proliferation experimental results are analyzed using Calcsyn 2.1 software to mathematically determine the IC ₅₀ values. The Wilcoxon Rank Sum Test is used to compare the IC ₅₀ values of Crenolanib for a given mutation.
Cell Research	Cell growth is measured by MTT. Briefly, cells were seeded at 2 × 10 ³ cells per well in 96-well culture plates in triplicate. After the addition of various concentrations of drugs, cells are incubated for 72 h. At the end of culture, the plates are washed twice with PBS,

A DRUG SCREENING EXPERT

Cell Research	and cells are incubated with 200 μ l of RPMI 1640 containing 10% FCS and 0.25 mg/ml of MTT at 37 °C for 3 h. The absorbance of each well is measured with Titer-Tech 96-well multiscanner at 570 nm. The viable cell number is proportional to the absorbance. (Only for Reference)
---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Solubility Information

Solubility	Ethanol: 55.6 mg/mL (100.1 mM), Sonication is recommended. DMSO: 45 mg/mL (81.01 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (3.6 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8003 mL	9.0014 mL	18.0028 mL
5 mM	0.3601 mL	1.8003 mL	3.6006 mL
10 mM	0.180 mL	0.9001 mL	1.8003 mL
50 mM	0.036 mL	0.180 mL	0.3601 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Na HJ, et al. J Immunol. 2004, 173(2), 1276-1283.
Kim HS, et al. Toxicol Appl Pharmacol. 2006, 215(3), 317-329.
Kim CK, et al. Mol Pharmacol. 2010, 78(1), 142-150

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481