


## G(8-O-4)S Glycerol

## Chemical Properties

CAS No. :	1313434-74-0
Formula:	C <sub>21</sub> H <sub>28</sub> O <sub>10</sub>
Molecular Weight:	440.44
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	G(8-O-4)S Glycerol is a natural product for research related to life sciences. The catalog number is TN5408 and the CAS number is 1313434-74-0.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2705 mL	11.3523 mL	22.7046 mL
5 mM	0.4541 mL	2.2705 mL	4.5409 mL
10 mM	0.227 mL	1.1352 mL	2.2705 mL
50 mM	0.0454 mL	0.227 mL	0.4541 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481