

Nevirapine dimer

Chemical Properties


CAS No. : 1391054-30-0

Formula: C₃₀H₂₆N₈O₂

Molecular Weight: 530.58

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nevirapine dimer is a non-nucleoside reverse transcriptase inhibitor (NNRTI).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8847 mL	9.4236 mL	18.8473 mL
5 mM	0.3769 mL	1.8847 mL	3.7695 mL
10 mM	0.1885 mL	0.9424 mL	1.8847 mL
50 mM	0.0377 mL	0.1885 mL	0.3769 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481