


cIAP1 Ligand-Linker Conjugates 1

Chemical Properties

CAS No. :	2357114-75-9
Formula:	C42H57N5O9S
Molecular Weight:	808.0
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	cIAP1 Ligand-Linker Conjugates 1 is composed of an IAP ligand for the E3 ubiquitin ligase and a PROTAC linker. This compound, cIAP1 Ligand-Linker Conjugates 1, is particularly useful in the development of SNIPERs[1].
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2376 mL	6.1881 mL	12.3762 mL
5 mM	0.2475 mL	1.2376 mL	2.4752 mL
10 mM	0.1238 mL	0.6188 mL	1.2376 mL
50 mM	0.0248 mL	0.1238 mL	0.2475 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Nello Mainolfi, et al. Irak degraders and uses thereof. US20190192668A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481