

Diocanoylglycol

Chemical Properties

CAS No. : 627-86-1


Formula: C₁₈H₃₄O₄

Molecular Weight: 314.46

Pure form: -20°C for 3 years | In solvent: -80°C for 1

Storage: year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Diocanoylglycol (ethane-1,2-diyl dioctanoate) is a diacylglycerol kinase inhibitor.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 18.33 mg/mL (58.29 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.36 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1801 mL	15.9003 mL	31.8005 mL
5 mM	0.636 mL	3.1801 mL	6.3601 mL
10 mM	0.318 mL	1.590 mL	3.1801 mL
50 mM	0.0636 mL	0.318 mL	0.636 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Slack BE, Richardson UI, Nitsch RM, Wurtman RJ. Dioctanoylglycerol stimulates accumulation of [methyl-14C] choline and its incorporation into acetylcholine and phosphatidylcholine in a human cholinergic neuroblastoma cell line. Brain Res. 1992 Jul 10;585(1-2):169-76.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481