


## Trigothysoid N

## Chemical Properties

CAS No. : 1501943-08-3  
 Formula: C44H58O13  
 Molecular Weight: 794.935  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Trigothysoid N may have anti-HIV activity.
Targets(IC50)	HIV Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.258 mL	6.2898 mL	12.5796 mL
5 mM	0.2516 mL	1.258 mL	2.5159 mL
10 mM	0.1258 mL	0.629 mL	1.258 mL
50 mM	0.0252 mL	0.1258 mL	0.2516 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481