

ATR-IN-11

Chemical Properties

CAS No. :

Formula: C₂₅H₃₀N₆O₂

Molecular Weight: 446.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ATR-IN-11 (Compound Hit01) is a potent inhibitor of ataxia capillaris and Rad3-related (ATR) kinases. As a promising lead compound for developing drugs targeting the ATR kinase, a key regulatory protein in the DNA damage response (DDR) that senses replication stress (RS), ATR-IN-11 has significant potential for cancer research.
Targets(IC50)	ATM/ATR,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2394 mL	11.1972 mL	22.3944 mL
5 mM	0.4479 mL	2.2394 mL	4.4789 mL
10 mM	0.2239 mL	1.1197 mL	2.2394 mL
50 mM	0.0448 mL	0.2239 mL	0.4479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481