


WE 1073

Chemical Properties

CAS No. : 88986-12-3
 Formula: C₁₅H₁₀BrClN₄O₂S
 Molecular Weight: 425.69
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	WE 1073 is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3491 mL	11.7456 mL	23.4913 mL
5 mM	0.4698 mL	2.3491 mL	4.6983 mL
10 mM	0.2349 mL	1.1746 mL	2.3491 mL
50 mM	0.047 mL	0.2349 mL	0.4698 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481