


Mabuterol hydrochloride

Chemical Properties

CAS No. :	54240-36-7
Formula:	C ₁₃ H ₁₉ Cl ₂ F ₃ N ₂ O
Molecular Weight:	347.20
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mabuterol hydrochloride is a selective β 2 adrenoreceptor agonist. Mabuterol has a specific effect on beta 2-adrenoceptors with no beta 1-stimulation.
Targets(IC50)	Others,Adrenergic Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8802 mL	14.4009 mL	28.8018 mL
5 mM	0.576 mL	2.8802 mL	5.7604 mL
10 mM	0.288 mL	1.4401 mL	2.8802 mL
50 mM	0.0576 mL	0.288 mL	0.576 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481