


Fmoc-Aib-OH

Chemical Properties

CAS No. :	94744-50-0
Formula:	C ₁₉ H ₁₉ N ₁ O ₄
Molecular Weight:	325.36
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Fmoc-Aib-OH (Fmoc-alpha-methylalanine) is an alanine acid derivative containing an Fmoc protecting group, which can be used for peptide synthesis and biopharmaceutical design.
Targets(IC50)	Amino Acids and Derivatives

Solubility Information

Solubility	DMSO: 180 mg/mL (553.23 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0735 mL	15.3676 mL	30.7352 mL
5 mM	0.6147 mL	3.0735 mL	6.147 mL
10 mM	0.3074 mL	1.5368 mL	3.0735 mL
50 mM	0.0615 mL	0.3074 mL	0.6147 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481