


Rwj 29009

Chemical Properties

CAS No. : 143164-10-7
 Formula: C₁₄H₁₈N₂O₅
 Molecular Weight: 326.37
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Rwj 29009 is a potassium channel activator with prominent coronary and peripheral vasodilating actions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.064 mL	15.320 mL	30.6401 mL
5 mM	0.6128 mL	3.064 mL	6.128 mL
10 mM	0.3064 mL	1.532 mL	3.064 mL
50 mM	0.0613 mL	0.3064 mL	0.6128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481