


Chlorindanol

Chemical Properties

CAS No. :	145-94-8
Formula:	C ₉ H ₉ ClO
Molecular Weight:	168.62
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Chlorindanol (Lanesta) is a new antiseptic agent. It is nonallergenic and appears to be suitable for use as a topical antiseptic or sanitizer.
Targets(IC50)	Others, Antibacterial

Solubility Information

Solubility	DMSO: 50 mg/mL (296.52 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (11.86 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9305 mL	29.6525 mL	59.3049 mL
5 mM	1.1861 mL	5.9305 mL	11.861 mL
10 mM	0.593 mL	2.9652 mL	5.9305 mL
50 mM	0.1186 mL	0.593 mL	1.1861 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

BAILEY JH, et al. J Am Pharm Assoc Am Pharm Assoc (Baltim). 1959 Apr;48(4):212-6.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481