


## cis-1,2-Dihydroxycyclohexane

## Chemical Properties

CAS No. : 1792-81-0  
 Formula: C<sub>6</sub>H<sub>12</sub>O<sub>2</sub>  
 Molecular Weight: 116.16  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	cis-1,2-Dihydroxycyclohexane is a diol and cyclohexane derivative with potential carbonic anhydrase inhibitory activity for biochemical experiments and drug synthesis research.
Targets(IC50)	Carbonic Anhydrase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.6088 mL	43.0441 mL	86.0882 mL
5 mM	1.7218 mL	8.6088 mL	17.2176 mL
10 mM	0.8609 mL	4.3044 mL	8.6088 mL
50 mM	0.1722 mL	0.8609 mL	1.7218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481