


## 4-Desmethyldecaline

## Chemical Properties

CAS No. :	54422-29-6
Formula:	C <sub>25</sub> H <sub>29</sub> NO <sub>5</sub>
Molecular Weight:	423.50
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-Desmethyldecaline is a biochemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3613 mL	11.8064 mL	23.6128 mL
5 mM	0.4723 mL	2.3613 mL	4.7226 mL
10 mM	0.2361 mL	1.1806 mL	2.3613 mL
50 mM	0.0472 mL	0.2361 mL	0.4723 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Hanaoka M, Ogawa N, Arata Y. Letter: Synthesis of (plus or minus)-desmethyldecaline. Chem Pharm Bull (Tokyo). 1974 Aug;22(8):1945-6. PubMed PMID: 4430050.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481