


3',5'-Bis-O-(triphenylmethyl)uridine

Chemical Properties

CAS No. : 4710-75-2
 Formula: C₄₇H₄₀N₂O₆
 Molecular Weight: 728.83
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	3',5'-Bis-O-(triphenylmethyl)uridine is a Nucleoside Derivative - Protected nucleoside with NH ₂ /OH group.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3721 mL	6.8603 mL	13.7206 mL
5 mM	0.2744 mL	1.3721 mL	2.7441 mL
10 mM	0.1372 mL	0.686 mL	1.3721 mL
50 mM	0.0274 mL	0.1372 mL	0.2744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481