

## 2-(2-hydroxypropan-2-yl)-4a,8-dimethyl-decahydronaphthalene-1,5-diol

## Chemical Properties

CAS No. :

Formula: C15H28O3

Molecular Weight: 256.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	2-(2-hydroxypropan-2-yl)-4a,8-dimethyl-decahydronaphthalene-1,5-diol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9003 mL	19.5015 mL	39.0031 mL
5 mM	0.7801 mL	3.9003 mL	7.8006 mL
10 mM	0.390 mL	1.9502 mL	3.9003 mL
50 mM	0.078 mL	0.390 mL	0.7801 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481