

Bayogenin methyl ester

Chemical Properties


CAS No. :

Formula: C₃₁H₅₀O₅

Molecular Weight: 502.736

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bayogenin methyl ester is a useful organic compound for research related to life sciences and the catalog number is T125298.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9891 mL	9.9455 mL	19.891 mL
5 mM	0.3978 mL	1.9891 mL	3.9782 mL
10 mM	0.1989 mL	0.9945 mL	1.9891 mL
50 mM	0.0398 mL	0.1989 mL	0.3978 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481