


Compound EN300-113927

Chemical Properties

CAS No. : 118049-54-0
 Formula: C₁₆H₁₁N₃O
 Molecular Weight: 261.278
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description Compound EN300-113927 is a high-purity biochemical reagent with CAS number 118049-54-0, suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8273 mL	19.1366 mL	38.2731 mL
5 mM	0.7655 mL	3.8273 mL	7.6546 mL
10 mM	0.3827 mL	1.9137 mL	3.8273 mL
50 mM	0.0765 mL	0.3827 mL	0.7655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481