


3-(1H-pyrrol-1-yl)aniline

Chemical Properties

CAS No. :	89353-42-4
Formula:	C10H10N2
Molecular Weight:	158.1998
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	3-(1H-pyrrol-1-yl)aniline, with CAS No. 89353-42-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3-(1H-pyrrol-1-yl)aniline provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.3211 mL	31.6056 mL	63.2111 mL
5 mM	1.2642 mL	6.3211 mL	12.6422 mL
10 mM	0.6321 mL	3.1606 mL	6.3211 mL
50 mM	0.1264 mL	0.6321 mL	1.2642 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481