


6-fluorinated-aristeromycin2c

Chemical Properties

CAS No. :	2236598-84-6
Formula:	C ₁₁ H ₁₃ F ₂ N ₅ O ₃
Molecular Weight:	301.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	6-fluorinated inated-aristeromycin2c showed antiviral activity against heat virus, Zika virus and Sars coronavirus.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3195 mL	16.5975 mL	33.195 mL
5 mM	0.6639 mL	3.3195 mL	6.639 mL
10 mM	0.332 mL	1.6598 mL	3.3195 mL
50 mM	0.0664 mL	0.332 mL	0.6639 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481