

Schineolignin B

Chemical Properties

CAS No. :	1352185-26-2
Formula:	C ₂₂ H ₃₀ O ₅
Molecular Weight:	374.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Schineolignin B is a useful organic compound for research related to life sciences and the catalog number is T126011.
Targets(IC50)	Others,HIV Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6704 mL	13.3522 mL	26.7044 mL
5 mM	0.5341 mL	2.6704 mL	5.3409 mL
10 mM	0.267 mL	1.3352 mL	2.6704 mL
50 mM	0.0534 mL	0.267 mL	0.5341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481