


4-Piperidinyl(1-pyrrolidinyl)Methanone hydrochloride

Chemical Properties

CAS No. : 188979-07-9
 Formula: C₁₀H₁₉ClN₂O
 Molecular Weight: 218.724
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.

HCl



Biological Description

Description	4-Piperidinyl(1-pyrrolidinyl)Methanone hydrochloride, with CAS No. 188979-07-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Piperidinyl(1-pyrrolidinyl)Methanone hydrochloride provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5721 mL	22.8603 mL	45.7206 mL
5 mM	0.9144 mL	4.5721 mL	9.1441 mL
10 mM	0.4572 mL	2.286 mL	4.5721 mL
50 mM	0.0914 mL	0.4572 mL	0.9144 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481