

Osthenone

Chemical Properties

CAS No. :	112789-90-9
Formula:	C14H12O4
Molecular Weight:	244.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Osthenone is a useful organic compound for research related to life sciences and the catalog number is T126372.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0943 mL	20.4717 mL	40.9433 mL
5 mM	0.8189 mL	4.0943 mL	8.1887 mL
10 mM	0.4094 mL	2.0472 mL	4.0943 mL
50 mM	0.0819 mL	0.4094 mL	0.8189 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481