


Taxumairol R

Chemical Properties

CAS No. :	244167-04-2
Formula:	C37H44O15
Molecular Weight:	728.74
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Taxumairol R is a natural product from Taxus mairei.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3722 mL	6.8612 mL	13.7223 mL
5 mM	0.2744 mL	1.3722 mL	2.7445 mL
10 mM	0.1372 mL	0.6861 mL	1.3722 mL
50 mM	0.0274 mL	0.1372 mL	0.2744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481