

Anti-MERS-2E6 mAb

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-MERS-2E6 mAb (MERS-2E6), a human IgG1 neutralizing antibody expressed in CHO cells, effectively competes with the virus Spike protein for binding to the CD26 receptor, thereby inhibiting viral invasion into host cells.
Targets(IC50)	SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481