

## Demethylmangostanin

## Chemical Properties


CAS No. :

Formula: C<sub>46</sub>H<sub>43</sub>O<sub>12</sub>

Molecular Weight: 787.838

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Demethylmangostanin is a useful organic compound for research related to life sciences and the catalog number is T125275.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2693 mL	6.3465 mL	12.6929 mL
5 mM	0.2539 mL	1.2693 mL	2.5386 mL
10 mM	0.1269 mL	0.6346 mL	1.2693 mL
50 mM	0.0254 mL	0.1269 mL	0.2539 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481