


WAG003

Chemical Properties

CAS No. : 1400302-65-9
 Formula: C₂₄H₂₈I₂NO₃
 Molecular Weight: 632.29
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAG003 is a Tudor domain inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5816 mL	7.9078 mL	15.8155 mL
5 mM	0.3163 mL	1.5816 mL	3.1631 mL
10 mM	0.1582 mL	0.7908 mL	1.5816 mL
50 mM	0.0316 mL	0.1582 mL	0.3163 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481