


## N-acetyl Dapsone (D4')

## Chemical Properties

CAS No. :	2070015-28-8
Formula:	C <sub>14</sub> H <sub>14</sub> N <sub>2</sub> O <sub>3</sub> S
Molecular Weight:	294.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-acetyl Dapsone-D4 is the deuterium-labeled metabolite of N-acetyl Dapsone (TMSM-1691).
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3972 mL	16.986 mL	33.972 mL
5 mM	0.6794 mL	3.3972 mL	6.7944 mL
10 mM	0.3397 mL	1.6986 mL	3.3972 mL
50 mM	0.0679 mL	0.3397 mL	0.6794 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481