


2-Methyl-6-(p-tolyl)heptane-2,3-diol

Chemical Properties

CAS No. :	117421-22-4
Formula:	C ₁₅ H ₂₄ O ₂
Molecular Weight:	236.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Methyl-6-(p-tolyl)heptane-2,3-diol is a natural product for research related to life sciences. The catalog number is TN2778 and the CAS number is 117421-22-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.231 mL	21.1551 mL	42.3101 mL
5 mM	0.8462 mL	4.231 mL	8.462 mL
10 mM	0.4231 mL	2.1155 mL	4.231 mL
50 mM	0.0846 mL	0.4231 mL	0.8462 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481