


Phenol, 3-iodo-

Chemical Properties

CAS No. :	626-02-8
Formula:	C ₆ H ₅ IO
Molecular Weight:	220.01
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenol, 3-iodo- is a bioactive chemical.
Targets(IC50)	Others,Endogenous Metabolite

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5452 mL	22.7262 mL	45.4525 mL
5 mM	0.909 mL	4.5452 mL	9.0905 mL
10 mM	0.4545 mL	2.2726 mL	4.5452 mL
50 mM	0.0909 mL	0.4545 mL	0.909 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481