


Dihydrolycorine

Chemical Properties

CAS No. :	6271-21-2
Formula:	C ₁₆ H ₁₉ NO ₄
Molecular Weight:	289.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	1. Dihydrolycorine-HCL(DL) shows hypotensive effects, it can block alpha 1-adrenoceptors. 2. Dihydrolycorine is an inhibitors of protein synthesis in eukarytic cells, it halts protein synthesis in eukaryotic cells by inhibiting the peptide bone formation step. 3. Dihydrolycorine and nimodipine protects against anoxia damage of brain in rat.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 2.9 mg/mL (10.02 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (3.46 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4563 mL	17.2813 mL	34.5626 mL
5 mM	0.6913 mL	3.4563 mL	6.9125 mL
10 mM	0.3456 mL	1.7281 mL	3.4563 mL
50 mM	0.0691 mL	0.3456 mL	0.6913 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Jimenez A , Santos A , Alonso G , et al. Inhibitors of protein synthesis in eukaryotic cells. Comparative effects of some amaryllidaceae alkaloids.[J]. BBA Section Nucleic Acids And Protein Synthesis, 1976, 425(3):342-348.

Shen L, Zhao J, Xia Y, et al. Lycorine derivative effectively inhibits the replication of coronaviruses both in vitro and in vivo. hLife. 2023

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481