

Compound 0080-0011

## Chemical Properties


CAS No. :

Formula: C<sub>21</sub>H<sub>24</sub>N<sub>2</sub>O<sub>2</sub>

Molecular Weight: 336.435

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0080-0011 is a useful organic compound for research related to life sciences and the catalog number is T131706.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9723 mL	14.8615 mL	29.723 mL
5 mM	0.5945 mL	2.9723 mL	5.9446 mL
10 mM	0.2972 mL	1.4861 mL	2.9723 mL
50 mM	0.0594 mL	0.2972 mL	0.5945 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481