

3,5,25-Trihydroxyergostan-6-one

Chemical Properties


CAS No. :

Formula: C₁₉H₃₀O₃

Molecular Weight: 306.446

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,5,25-Trihydroxyergostan-6-one is a useful organic compound for research related to life sciences and the catalog number is T131391.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2632 mL	16.3159 mL	32.6318 mL
5 mM	0.6526 mL	3.2632 mL	6.5264 mL
10 mM	0.3263 mL	1.6316 mL	3.2632 mL
50 mM	0.0653 mL	0.3263 mL	0.6526 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481