


Shiraiachrome A

Chemical Properties

CAS No. :	124709-39-3
Formula:	C ₃₀ H ₂₆ O ₁₀
Molecular Weight:	546.52
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Shiraiachrome A is a natural product for research related to life sciences. The catalog number is T5525 and the CAS number is 124709-39-3.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8298 mL	9.1488 mL	18.2976 mL
5 mM	0.366 mL	1.8298 mL	3.6595 mL
10 mM	0.183 mL	0.9149 mL	1.8298 mL
50 mM	0.0366 mL	0.183 mL	0.366 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481