


S101

## Chemical Properties

CAS No. : 362508-67-6  
 Formula: C<sub>24</sub>H<sub>23</sub>N<sub>3</sub>O<sub>3</sub>  
 Molecular Weight: 401.46  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	S101 is an inhibitor of proliferating T-cells, rescuing mice from superantigen-induced shock.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4909 mL	12.4545 mL	24.9091 mL
5 mM	0.4982 mL	2.4909 mL	4.9818 mL
10 mM	0.2491 mL	1.2455 mL	2.4909 mL
50 mM	0.0498 mL	0.2491 mL	0.4982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481