


Beta-Amyrenonol acetate

Chemical Properties

CAS No. :	5356-56-9
Formula:	C32H50O3
Molecular Weight:	482.74
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Beta-Amyrenonol acetate is a natural product for research related to life sciences. The catalog number is TN3408 and the CAS number is 5356-56-9.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0715 mL	10.3575 mL	20.7151 mL
5 mM	0.4143 mL	2.0715 mL	4.143 mL
10 mM	0.2072 mL	1.0358 mL	2.0715 mL
50 mM	0.0414 mL	0.2072 mL	0.4143 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481