


2-Phenylglycine

Chemical Properties

CAS No. :	2835-06-5
Formula:	C ₈ H ₉ NO ₂
Molecular Weight:	151.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Phenylglycine (DL-α-Phenylglycine) is a metabolite present in breast milk during lactation.
Targets(IC50)	Endogenous Metabolite

Solubility Information

Solubility	H ₂ O: 3.33 mg/mL (22.03 mM), when pH is adjusted to 12 with NaOH. Sonication is recommended. DMSO: < 1 mg/mL (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6155 mL	33.0775 mL	66.1551 mL
5 mM	1.3231 mL	6.6155 mL	13.231 mL
10 mM	0.6616 mL	3.3078 mL	6.6155 mL
50 mM	0.1323 mL	0.6616 mL	1.3231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Alexandre-Gouabau MC, et al. Comprehensive Preterm Breast Milk Metabotype Associated with Optimal Infant Early Growth Pattern. *Nutrients*. 2019 Feb 28;11(3).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481