


PFI-4

Chemical Properties

CAS No. : 900305-37-5
 Formula: C₂₁H₂₄N₄O₃
 Molecular Weight: 380.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	PFI-4 is a potent and selective and cell-permeable BRPF1 bromodomain inhibitor.
Targets(IC50)	Epigenetic Reader Domain
In vitro	PFI-4 is cell permeable, and acts by displacing the BRPF1 bromodomain from Histone H3.3. [1]

Solubility Information

Solubility	H ₂ O: < 1 mg/mL (insoluble or slightly soluble), DMSO: 14 mg/mL (36.8 mM),Sonication is recommended. Ethanol: 5 mg/mL (13.14 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (2.63 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6285 mL	13.1427 mL	26.2854 mL
5 mM	0.5257 mL	2.6285 mL	5.2571 mL
10 mM	0.2629 mL	1.3143 mL	2.6285 mL
50 mM	0.0526 mL	0.2629 mL	0.5257 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Demont EH, et al. ACS Med Chem Lett. 2014, 5(11), 1190-1195.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481