

Somatostatin 1-28 acetate

Chemical Properties

CAS No. :

Formula:


Molecular Weight:

Storage:

Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Somatostatin 1-28 acetate belongs to natural products, serving as a multi-target agonist with endocrine inhibitory activity. Somatostatin 1-28 acetate circulates in human plasma, accounting for approximately 46% of total somatostatin-like immunoreactivity, and can be used for studying hormone secretion regulatory mechanisms.
Targets(IC50)	Somatostatin

Solubility Information

Solubility	DMSO: 250 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Reference

Baldissera FG, Munoz-Perez MA, Holst JJ. Somatostatin 1-28 circulates in human plasma. Regul Pept. 1983 Apr;6(1):63-9.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481