

## GABAA receptor agent 6

## Chemical Properties


CAS No. : 1808463-81-1

Formula: C<sub>18</sub>H<sub>25</sub>N<sub>3</sub>O<sub>2</sub>

Molecular Weight: 315.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GABAA receptor agent 6 (compound 2027) is a potent $\gamma$ -GABAAR antagonist with low cellular membrane permeability ( $K_i = 0.56 \mu\text{M}$ ) [1].
Targets(IC50)	Others,GABA Receptor
In vitro	Compound 2027, a GABAA receptor agent, exhibits antagonist activity, with IC <sub>50</sub> values of 6.68 $\mu\text{M}$ for $\alpha 1\beta 2\delta$ , 1.03 $\mu\text{M}$ for $\alpha 4\beta 1\delta$ , 0.36 $\mu\text{M}$ for $\alpha 4\beta 2\delta$ , 4.13 $\mu\text{M}$ for $\alpha 6\beta 2\delta$ , 4.96 $\mu\text{M}$ for $\alpha 1\beta 2\gamma 2$ , 2.96 $\mu\text{M}$ for $\alpha 2\beta 2\gamma 2$ , 0.29 $\mu\text{M}$ for $\alpha 3\beta 2\gamma 2$ , and 0.59 $\mu\text{M}$ for $\alpha 5\beta 2\gamma 2$ subtypes, respectively [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1705 mL	15.8524 mL	31.7048 mL
5 mM	0.6341 mL	3.1705 mL	6.341 mL
10 mM	0.317 mL	1.5852 mL	3.1705 mL
50 mM	0.0634 mL	0.317 mL	0.6341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481