

3-Nitro-L-tyrosine ethyl ester hydrochloride

Chemical Properties


CAS No. :

Formula: C₁₁H₁₅ClN₂O₅

Molecular Weight: 290.7

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Nitro-L-tyrosine ethyl ester hydrochloride [catalog number T124433] is a valuable compound for life sciences research.
Targets(IC50)	Others,Amino Acids and Derivatives

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.440 mL	17.1999 mL	34.3997 mL
5 mM	0.688 mL	3.440 mL	6.8799 mL
10 mM	0.344 mL	1.720 mL	3.440 mL
50 mM	0.0688 mL	0.344 mL	0.688 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481