

CH₂COOH-PEG₃-CH₂COOH

Chemical Properties

CAS No. : 32775-08-9


Formula: C₁₀H₁₈O₈

Molecular Weight: 266.25

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CH ₂ COOH-PEG ₃ -CH ₂ COOH is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, consisting of two distinct ligands linked together—one for an E3 ubiquitin ligase and the other for a target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7559 mL	18.7793 mL	37.5587 mL
5 mM	0.7512 mL	3.7559 mL	7.5117 mL
10 mM	0.3756 mL	1.8779 mL	3.7559 mL
50 mM	0.0751 mL	0.3756 mL	0.7512 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481