

N1-(N,N-Dimethylaminocarbonyl)-pseudouridine

Chemical Properties

CAS No. :


Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N1-(N,N-Dimethylaminocarbonyl)-pseudouridine is a nucleoside derivative classified as a C-nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481