


Celosin J

Chemical Properties

CAS No. :	1623405-29-7
Formula:	C ₅₈ H ₉₀ O ₂₈
Molecular Weight:	1235.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Celosin H, celosin I, celosin J could be used as chemical markers for the quality control of <i>C. argentea</i> seeds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8095 mL	4.0475 mL	8.0951 mL
5 mM	0.1619 mL	0.8095 mL	1.619 mL
10 mM	0.081 mL	0.4048 mL	0.8095 mL
50 mM	0.0162 mL	0.081 mL	0.1619 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481