

VGX100

Chemical Properties

CAS No. :

Formula:

Molecular Weight: 143.58 kDa

Store at low temperature

Storage: -20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | VGX100 is a humanized antibody targeting VEGFC for the study of solid tumors. |
| Targets(IC50) | VEGFR |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 0.007 mL | 0.0348 mL | 0.0696 mL |
| 5 mM | 0.0014 mL | 0.007 mL | 0.0139 mL |
| 10 mM | 0.0007 mL | 0.0035 mL | 0.007 mL |
| 50 mM | 0.0001 mL | 0.0007 mL | 0.0014 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481