

RO5458640

Chemical Properties

CAS No. :

Formula:

Molecular Weight: 146.02 kDa

Store at low temperature

Storage: -20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	RO5458640 (RG-7212) is a humanized antibody targeting TWEAK, used in tumor research.
Targets(IC50)	Others
In vitro	RO5458640 is a human antibody targeting TNFSF12/TWEAK that is involved in the regulation of inflammatory processes in neurons and tumor tissues. [1]

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.0342 mL	0.0685 mL
5 mM	0.0014 mL	0.0068 mL	0.0137 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481