

## Ladiratuzumab

## Chemical Properties

CAS No. :	1629760-28-6
Formula:	
Molecular Weight:	146.1 kDa
Storage:	Store at low temperature -20°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	LadiratuzumAb (hLIV22) is a humanized monoclonal antibody targeting LIV-1/ZIP6. LadiratuzumAb is a naked antibody to SGN-LIV1A and is capable of synthesizing an antibody-activated molecule coupling (ADC), SGN-LIV1A.LadiratuzumAb can be used for the prevention and treatment of breast cancer.
Targets(IC50)	Others,ADC Antibody

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.0342 mL	0.0684 mL
5 mM	0.0014 mL	0.0068 mL	0.0137 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481