

Eltivutabart

Chemical Properties

CAS No. : 2935496-38-9

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Eltivutabart is a humanized monoclonal antibody targeting ENTPD1/CD39. It increases extracellular ATP and reduces immunosuppressive adenosine, thereby promoting antitumor immunity. It is suitable for studying metastatic pancreatic cancer.
Targets(IC50)	Phosphatase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481