

GSK3050002

Chemical Properties

CAS No. :

Formula:

Molecular Weight: 144.24 kDa

Store at low temperature

Storage: -20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	GSK3050002 (E-6071) is a humanized antibody targeting CCL20, which can be used to study immune system diseases and digestive system diseases.
Targets(IC50)	Others
In vivo	Aggregation of GSK3050002 was time-dependent, and cynomophores receiving 30 and 300 mg/kg/ week of GSK3050002 subcutaneously or intravenous showed increased dose-responsiveness in the incidence and severity of inflammation in multiple organs in the monkeys. [1]

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0069 mL	0.0347 mL	0.0693 mL
5 mM	0.0014 mL	0.0069 mL	0.0139 mL
10 mM	0.0007 mL	0.0035 mL	0.0069 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481