

Tanrurubart

Chemical Properties

CAS No. : 2065212-40-8

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tanrurubart is a humanized monoclonal antibody targeting C1q in the classical complement pathway, applicable for studying Guillain-Barré syndrome (GBS) by preventing persistent neuroinflammation and nerve damage in the early stages of GBS.
Targets(IC50)	Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481